


# Tap Adapters

The T-safe Tap Adapter system enables secure connection of the T-safe Medical Tap Filter to tap outlets via a push fit mechanism. A wide range of Tap Adapters are available in variable configurations and thread diameters, with a Universal Connector for non-threaded tap outlets.



## Order Information

Tap Adapter	Code	Item	Code
Tap Filter Adapter M24m / M22f, CV, 6LFR	90-001804	Adapter Connector 21.5mm m to M24f	01-990198
Tap Filter Adapter M24m / M22f, Long, CV, 6LFR, FS	90-001806	Adapter Connector ½" f to M24f	90-001009
Tap Filter Adapter M24m, O-Ring, CV, 6LFR	90-001814	Adapter Connector ½" m to M24f	01-990248
Tap Filter Adapter M24m, Long, O-Ring, CV, 6LFR, FS	90-001816	Adapter Connector M18m to M24f	01-990247
		Universal Adapter Connector 15-29mm	90-001851

Key: (m) male (f) female CV: Check Valve FR: Flow Regulator FS: Flow Stop